

SMV025

Type IIR Face Mask with Visor



The surgical Face Mask with Visor is recommended for use wherever there is risk of fluid contamination, for example during splash or spray prone procedures. Constructed from four layers to differentiate them from standard masks, these Self Protection masks offer fluid resistance of 120mm/Hg. The masks fasten securely with ties at the nape and crown of your head, and offer an anti-fog visor that provides an enhanced level of protection when compared with conventional IIR face masks.

Key Features

- › Splash Resistant
- › Comfortable fit
- › Latex free
- › High Performance 4 ply filter
- › Anti-fog Visor
- › Tie back fastening
- › Non-sterile

Quality

Our Face Mask with Visor is manufactured under EN 14683 IIR and ASTM F2100 LEVEL 3. This ensures that we deliver products and services that consistently meet customer requirements, guaranteeing full customer satisfaction.

**ASTM
F2100
LEVEL 3**




Warning

- › This product cannot prevent inhalation of harmful or poisonous gases
- › This product is disposable and cannot be reused


Material Composition

Outer 2 Layers:	Spunbond Polypropylene
Filter:	Meltblown Polypropylene
Inner Layer:	Spunbond Polypropylene
Nose Bridge:	Aluminium Wire
Visor:	PVC
Ties:	Polypropylene




4 Layer Non-woven Fabrics


Instructions for Use




1 Make a dent in the middle of the nose clip and curve the visor.



2 Fasten the upper ties at the crown of your head.



3 Extend the mask to fully cover your mouth, nose, and chin. Mold the nose clip over your nose bridge.



4 Fasten the lower ties at the nape.

Sizing & Packaging Details

- › The Self Protection Surgical Mask with Visor must be stored in a clean, dry environment and avoid direct sunlight.
- › Suggested shelf life: 24 months

✓ 25 pcs/box

✓ 600 pcs/carton

✓ Dimensions: 175 x 302 mm(Mask: 175 x 95 mm)

Toffeln®

toffeln.com

Unit 2 Bridge Road, Kingswood, Bristol,
BS15 4FW, UK

✉ sales@toffeln.com

☎ 0345 500 4433