

TPPCW100

Community Wound Care Procedure Pack

The TPPCW100 Community Wound Care Procedure Pack is designed to provide all the essential tools and accessories needed for wound care. This pack includes 7 high-quality components, each tailored for efficient and safe community wound care.

See Pack Components for a bill of materials.



Features

- » **Comprehensive Kit:** Each procedure pack includes 7 components tailored to treat wounds.
- » **Sterilised and Ready for Use:** All components are sterilised and double-packed in a sterile environment ready for immediate use.
- » **High Quality:** Made from durable, medical-grade materials to ensure lasting performance and reliability in the operating theatre.
- » **Compliance:** Meets all Health & Safety regulations and is manufactured to ISO 13485:2016 standards.

Benefits

- » **Enhanced Efficiency:** Having all the necessary tools in one pack streamlines preparation, saving your team time and boosting productivity.
- » **Patient Safety:** A pre-sterilised procedure pack reduces the risk of patient infection.
- » **More Storage Space:** Using one procedure pack significantly improves hospital storage space.
- » **Reduce Waste:** Consolidating all components significantly reduces packaging waste.
- » **Cost Effective:** Combining essential tools into one pack is more economical than purchasing items individually.



Supply Chain
Product

Pack Components

Product	Size	Qty
Apron	Standard	1
Dressing towel	Standard	1
Swabs non-woven	7.5 x 7.5cm	3
Waste bag	Standard	1
Wrapping Drape	50 x 50cm	1



Instructions of Use

1. Inspect the outer packaging for damage or tampering, check the sterility indicator, and verify the expiration date before opening.
2. Thoroughly wash and sanitise your hands according to standard medical protocols.
3. Open the sterile pack in a clean and sterile environment, following the medical protocol to avoid contamination.
4. Carefully remove the tools and accessories from the pack and arrange them on a sterile surface in the order of use.
5. Ensure all components listed in the pack contents are present and in good condition before beginning the procedure.
6. Use each tool as required by the specific medical procedure, adhering to standard operating protocols and guidelines.
7. Safely dispose items after the procedure.
8. Document the use of the procedure pack and any observations or anomalies during the procedure for future reference.

Packaging & Storage

Dimensions per: Pack: W 14 x L 16 x H 2 (cm)
Carton: W: 40 x L 60 x H 32 (cm)

Weight per: Pack: 50g
Carton: 3kg

Storage: Store in a dry environment away from direct sunlight. If stored correctly, each procedure pack has a 5-year shelf life from the date of manufacture.

Ordering Information

Product Code	NHS NPC	Qty (per carton)
TPPCW100	EVS87298	60



Supply Chain Product